

Week 1: Introduction to BIOPAC

BIOE 320 Systems Physiology Laboratory

Objectives

1. To become familiar with the format of data display in the BIOPAC Student Lab data window
2. To learn how to position data within the window by using software tools and pull down menus
3. To learn how to select and use correct measurement tools for extracting information from the data
4. To learn how to use journal to record measurements and write notes

Background

The lab is based on the BIOPAC Student Lab System that allows data acquisition and analysis of a variety of physiological signals. **Hardware** components include electrodes, input transducers, connecting cables, and the acquisition unit (MP30, MP35, or MP36) shown in Fig. 1. The **software** used in this lab consists of the BIOPAC Student Lab Lessons and the BIOPAC Student Lab PRO pre-installed on your computer.




Figure 1: Schematic of MP3x acquisition unit and attached transducer

The lab will introduce you to the BIOPAC system. To complete the in-lab activities, read the instructions carefully and perform each of the steps described in the exercises while answering the questions listed below.



Info: Your knowledge and comfort level in manipulating the data obtained using BIOPAC are critical for success in future labs.

Getting Started

Before beginning the lab activity, open the BIOPAC tutorial that can be found in: (need to figure out links)

1. Turn on the MP3X acquisition unit by flipping the switch found on the back panel.




Figure 2: MP3X acquisition unit - front panel (top) and back panel (bottom) are shown

2. Launch the BIOPAC Student Lab Lessons Software (*Start* → *All Programs* → *BIOPAC Student Lab Lessons 3.7.6*)

⚠ **Warning:** If you receive the warning prompt shown in Fig. 3, make sure the acquisition box is on, check the connections to the computer, and restart the hardware and software. If the problem persists, call the instructor or a TA.




Figure 3: Connection error prompt

3. Scroll down the lesson menu to "Review Saved Data," select it, and click OK. A window resembling Fig. 4 (left) should appear.




Figure 4: Opening data files

4. If the default folder does not come up, navigate to 'C:\Program Files\BIOPAC Systems, Inc.\BIOPAC Student Lab 3.7\BSL Lessons 3.7.6\Data Files\Sample Data\SampleData-L07'
5. Select and open the "Sample Data" folder. A window resembling Fig. 4 (right) should appear.
6. Open "SampleData-L07".
7. The monitor screen should display a data window showing ECG (top waveform) and Pulse (bottom waveform) data (Fig. 5)




Figure 5: L07 data file

Data Viewing Tools

1. The software consists of 2 main sections: data window (area with graphs) and journal window (area with text notes). Make sure you can find these on your computer.
2. In this sample data, the data window contains two series of simultaneously recorded waveforms: an electrocardiogram (ECG) and a pulse. These waveforms depend on

the physiological parameter measured and will change depending on the lesson recorded.

3. Select the different markers (found above the marker bar in blue diamonds) and note from the markers what experiments were performed with this data set.
4. Turn on and off channel 1 (ECG) and channel 40 (pulse). Remember that the channel boxes (found above the marker bar on the left) allow you to display or hide channels from view in order to concentrate or print out only the desired waveforms. The disabled channel will be hidden from view, but can be turned back on at any time.

Data Analysis Tools: Selection, I-beam, and Zoom Tools ()

1. Use the selection tool (arrow) to better select the section from 3 to 40 seconds: left click on the x-axis (seconds) to open the prompt seen in Fig. 6 and change the scale range to display the desired data section.




Figure 6: Changing the display for the horizontal scale

2. Use the I-beam tool (looks like an I) to highlight different sections of the data. Notice that the measurement boxes (Delta T, BPM, P-P) will always display information corresponding to the highlighted region.
3. Use the zoom tool (magnifying glass) to select an area of the data to be enlarged so as to make some analyses easier to perform. Use the zoom tool to display the waveform at approximately 12.6 seconds. Hold down the mouse button and drag the cursor downward and to the right so that the selected area includes approximately six spikes (see Fig. 7).




Figure 7: Selecting 6 cycles using the zoom tool

4. Release the mouse button. Your data should not look like Fig. 8. If you selected the wrong part of the data to magnify, pull down the Display menu and choose "Zoom Previous", then try again.




Figure 8: Selected region using zoom tool

Channel Measurement Box

1. Right above the ECG recording channel, you will find the measurement boxes (Delta T, BPM, P-P). The first box indicates the *recording channel* selected, the second box indicates the *function or measurement* to be taken, and the third box contains the *measurement value* (see Fig. 9).




Figure 9: Channel measurement boxes

2. Using the first cycle (first beat), answer the following questions:
 - (a) What is your measured value for delta T?
 - (b) What does delta T represent?
 - (c) What is your measured value for BPM?
 - (d) What does BPM represent?
 - (e) What are the "max" and "min" for the waveform on ECG?

- (f) What are the "max" and "min" for the waveform on Pulse?
- (g) What do the "max" and "min" values represent?
- (h) Try out the options for "area", "P-P", "integral", "stdev", and two additional ones (see Table [add table at end] for available options and explanations). Explain in your own words what they are and the resulting value.

Overlap/Split

- These buttons allow you to toggle the display style between an oscilloscope and a chart recorder (see Fig. 10)




Figure 10: Display types: split (top) and overlap (bottom)



- Click the overlap button to activate the type of display. Once the signals are overlapped, try to move each signal up and down separately for alignment.

i **Info:** Select the channel box to select the appropriate channel to adjust.

- Return the display type to split.

Markers

Markers are used to reference important events (and their associated locations) in the data. There are two types of markers:

-  *Append markers*: These appear as a diamond above the marker text box and are blue when active. Append markers are automatically inserted when you begin each new recording segment and are marked with time data. Markers can also be manually entered during recording by pressing F9.
-  *Event markers*: These appear as inverted triangles below the marker text region and are yellow when active. You may add markers to your data after it has been recorded simply by clicking within the marker region using the selection tool. This new marker will then become the current active marker, and you may type in the marker text box.




Figure 11: Marker tools are located in the top right, indicated by the arrows

Using the marker tools on the right edge of the marker region (Fig. 11), you may:

- Change the active marker using the left and right pointing arrows
- Choose from the pop-up menu to find, clear, show, or import to journal certain or all markers (Fig. 12).




Figure 12: Marker tools pop-up menu

1. How do you insert a marker while data are being recorded?
2. How do you insert text for a marker?
3. How do you insert an event marker after data have been recorded? Do this as practice.